

QUICK TIPS FOR THE CORRECT USE OF THE SAFETY ARM SYSTEM™



LOWERING THE SAFETY ARM

- PULL AND HOLD THE SAFETY PIN
- LIFT THE ARM IN ITS ANGLE
- LOWER THE ARM INTO ITS HORIZONTAL POSITION
- RELEASE THE SAFETY PIN

RAISING THE SAFETY ARM

- PULL AND HOLD THE SAFETY PIN
- RAISE THE ARM INTO VERTICAL POSITION
- LOWER THE ARM INTO ITS ORIGINAL POSITION
- RELEASE THE SAFETY PIN

REMOVING THE SAFETY ARM

- WITH THE ARM IN VERTICAL POSITION, REMOVE THE LOCK PIN
- PULL AND HOLD THE SAFETY PIN
- REMOVE THE ARM
- RELEASE THE SAFETY PIN



Always refer to the product user guide for the correct installation and use of the equipment. Find it here:

© 2021, Technimount EMS Inc. — All Rights Reserved.

SAFETY ARM SYSTEM™ CONFIGURATION OPTIONS

INCREASED FLEXIBILITY FOR SPACE MAXIMIZATION AND

MULTIPLE DEVICE USE FOR 911 AND CRITICAL CARE TRANSPORT



TECHNIMOUNT
EMS™

Safety Arm Right Side
and IV Pro Adapter™



Safety Arm Right Side
and Safety Arm Angular Left Side



Safety Arm Angular



The **Safety Arm System™** is a secure stretcher mounting system designed to transport various types of medical devices, such as monitors, defibrillators and ventilators, which enables equipment to follow the patient during care.

With multiple configuration options available on the left and right side of the stretcher, the **Safety Arm System™** improves medical device management while increasing safety for EMS providers and patients during EMS and Critical Care Transport.



www.penncare.net
sales@penncare.net
1-800-392-7233