

## Product Comparison

LIFEPAK® 1000  
Defibrillator



LIFEPAK CR® Plus  
Defibrillator



USE		
Device Users	EMS, Fire, Police, Industrial or Frequent Users (use AED more than once a year)	Workplace or Infrequent Users (use AED less than once a year)
THERAPY FEATURES		
Defibrillation	Escalating Energy 200-300-360J	Escalating Energy 200-300-360J
Shock	Semi-automatic	Fully automatic
Manual Shock Option	Yes	No
DEVICE FEATURES		
Voice Prompts	Yes	Yes
Readiness Display	OK, BATTERY, SERVICE	OK, BATTERY, ATTENTION, SERVICE
Electrode Options	EDGE System™ Electrodes with REDI-PAK™ Preconnect System. Infant/Child Reduced Energy Defibrillation Electrodes.	QUIK-PAK™ Electrodes with REDI-PAK Preconnect System. Infant/Child Reduced Energy Defibrillation Electrodes.
Dust/Water Rating	IP55 (Protection from jetting water)	IPX4 (Protection from splashing water)
Post-Event Review	Infrared wireless to CODE-STAT™ software	Infrared wireless to CODE-STAT software
Power	Battery - 17 hours monitoring or 440 shocks	SafeGuard™ Power System - internal battery with replaceable battery charger - 210 minutes or 30 shocks
Included Accessories	Carry Case, Shoulder Strap, Two sets of electrodes, One battery	Carry Case, Two sets of electrodes, CHARGE-PAK™ battery charger, Ambu® Res-Cue Mask® Kit
Warranty	5 years	8 years

LIFEPAK defibrillators require a prescription in the U.S. Please consult your physician. Users should be trained in CPR and use of the defibrillator.

---

DISTRIBUTED BY:



©2012 Physio-Control, Inc. GDR 3311286\_A

 Physio-Control, Inc., 11811 Willows Road NE, Redmond, WA 98052 USA